

Post-it® Poster Paper (pages 16–17)

Assortability:

Not Available.

Colors for Standard Items:

White card stock with removable liner.

Copy Change:

- **\$25.00 (G)** per copy change—catalog quantities only.
- **Digital:** Copy change qualifies for any order of more than one design (artwork) grouped together. Minimum copy change quantity is 50. Minimum order quantity is 250.
- **Offset:** Two-color printing changes in copy allowed with only one of the two colors selected. Copy change location and color must remain the same. Minimum copy change quantity is 250. Minimum order quantity is 500.

Imprint Area:

Visit www.3M.com/templates for printing specifications.

Ink Charges:

No additional charge for PMS Inks except fluorescent, special mix and non-standard metallic **\$285.00 (G)** per color. Exact match to PMS ink color cannot be guaranteed. PMS Inks not available for digital printing.

Ink Colors:

- **Standard Ink Colors for Exact Color Match printing, no additional charge:**

Black	Yellow/Process Yellow
Green/PMS 347	Purple/PMS 527
Red/PMS 185	Blue/PMS 300
Brown/PMS 168	Teal/PMS 320
	Gold Metallic/PMS 872
	Silver Metallic/PMS 877
	Burgundy/PMS 215

NOTE: PMS number listed next to Standard Color is a very close representation of standard color. Exact match cannot be guaranteed. Not all colors provide high contrast on product. PMS, fluorescent, special mix and metallic inks are not available for digital printing.

Minimum Order:

- **11" x 17" Poster Paper (PP1117):** 250
 - **8-1/2" x 11" Poster Paper (PP811):** 100
 - **UNPRINTED Poster Paper (all sizes):** 50
- Call for quote.

Packaging Charges and Requirements:

- **Case quantity/weights:** See individual product pages.
- **Special shrink-wrapping or packaging:** Available at extra charge — call for quote.
- **Special packaging requirements:** Any special packaging requirements must accompany the order to ensure that an accurate delivery schedule can be established. Inserts or labels must arrive at the factory within **three (3)** working days of a complete order/proof approval, or the delivery schedule will be affected.

Plates, Screens and Set-up Charges:

No charge.

Printers for Blank Stock:

Poster Paper is designed for digital presses that can accommodate 270 GSM stock. **Poster Paper is not intended for use on inkjet or desktop printers.**

Production Time:

Normal production time is **five (5)** working days after receipt of complete order or after approval of proof. Call for production time on quantities larger than 4999.

Proof Charges:

Available on request with order. When a proof is requested, additional production time may be necessary. Call for confirmation of final delivery schedule.

- **E-Proof: \$8.00 (G).** Allow **three (3)** working days.
- **Black and White Faxed Proof: \$8.00 (G).**

Press Production:

- **Digital Press Proof: \$100.00 (G).**
 - **Offset Press Proof: Call for quote.**
- To assure customer satisfaction, we recommend a product proof for all 4-color process orders to ensure color expectations. This process takes **five (5)** working days to produce. Otherwise, provide a high quality color print of which is an accurate representation of the color expectations.

Prototypes:

Refer to Speculative Samples on page 40.

Restrictions:

- No bleeds.
- No imprint on liner.
- White paper only.
- No PMS, fluorescent, special mix or metallic inks for digital printing. All colors will be converted for digital printing. ■

Post-it® Pop-Up Note and Flag Dispensers (pages 28–29)

Assortability:

Not Available.

Colors for Standard Items:

- **Pop-Up Note Desktop Dispenser (DS330):**
 - Black base with Tinted Clear top
 - White paper insert.
 - White card stock lenticular insert.
 - White or Canary Yellow paper available only for printed notes.
 - Canary Yellow paper available only for unprinted notes.
- **Pop-Up Note and Flag Dispenser (DS100):**
 - Black base with Tinted Clear top
 - White or Canary Yellow paper available only for printed notes.
 - Canary Yellow paper available only for unprinted notes.
 - Red and Blue flags.

Copy Change:

Not Available.

Imprint Area:

Visit www.3M.com/templates for printing specifications.

Ink Charges:

No additional charge for PMS Inks except fluorescent, special mix and non-standard metallic **\$285.00 (G)** per color for all quantities. Exact match cannot be guaranteed.

Ink Colors:

- **Standard Ink Colors for Dispensers (DS330, DS100):**

Black	Yellow/Process Yellow
Purple/PMS 527	Green/PMS 347
Blue/PMS 300	Gold Metallic/PMS 872
Red/PMS 185	Silver Metallic/PMS 877
Brown/PMS 168	Burgundy/PMS 215
Teal/PMS 320	White

NOTE: PMS number listed next to Standard Color is a very close representation of standard color. Exact match cannot be guaranteed. Not all colors provide high contrast on products.

- **Paper and Lenticular Insert for DS330:** Printed digitally. All PMS inks will be converted for digital printing. Exact match to PMS ink color *cannot* be guaranteed. No fluorescent, special mix or metallic inks. **NOTE:** Colors of paper insert will shift due to tinted top.

Laser Engraving:

Laser engraving will be done on the underside of the clear plastic top. Listed catalog prices include **one (1)** engraving in **one (1)** location on **one (1)** side of note opening. For more than one laser engraving or laser engraving on both sides of note opening, call for quote and production time.

Minimum Order:

- **Pop-Up Note Desktop Dispenser (DS330):** 12 units
- **Pop-Up Note Desktop Dispenser with Lenticular Insert Card (DS330):** 1000 units
- **Pop-Up Note and Flag Dispenser (DS100):** 10 units

Packaging Charges and Requirements:

- **Case quantity/weights:** See individual product pages.
- **Pop-Up Note Desktop Dispenser (DS330):**
 - Packaged with one 3" x 3" 50-sheet, Canary Yellow, unprinted pop-up note pad.
- **Pop-Up Note and Flag Dispenser (DS100):**
 - Packaged with one 3" x 3" 50-sheet, Canary Yellow, unprinted pop-up note pad and one Red 1" flag dispenser and one Blue 1" flag dispenser.
- **Special shrink-wrapping or packaging:** Available at extra charge — call for quote.
- **Special packaging requirements:** Any special packaging requirements must accompany the order to ensure that an accurate delivery schedule can be established. Inserts or labels must arrive at the factory within **three (3)** working days of a complete order/proof approval, or the delivery schedule will be affected.

Paper Colors:

- **Pop-Up Note Desktop Dispenser (DS330):**
 - Unprinted: Canary Yellow, 50-sheet pop-up note pad.
 - Printed: Canary Yellow or White, 50-sheet pop-up note pad.
- **Pop-Up Note and Flag Dispenser (DS100):**
 - Unprinted: Canary Yellow, 50-sheet pop-up note pad.
 - Printed: Canary Yellow or White, 50-sheet pop-up note pad.

Plates, Screens and Set-up Charges:

No charge. **Note:** Two-sided sheet printing is not available.

Production Time:

On orders of fewer than 1499 dispensers, normal production time is **eight (8)** working days after receipt of complete order or after approval of proof. Add **five (5)** working days for sheet printing. On orders of fewer than 1499 dispensers, normal production time is **five (5)** working days for undecorated dispenser. Call for production time on lenticular insert cards for DS330.

Proof Charges:

Available on request with order. When a proof is requested, additional production time may be necessary. Call for confirmation of final delivery schedule.

- **E-Proof: \$8.00 (G).** Allow **three (3)** working days.
- **Black and White Faxed Proof: \$8.00 (G).**

Press Production:

- To assure customer satisfaction, we recommend a product proof for all 4-color process orders to ensure color expectations. This process takes **five (5)** working days to produce — for pricing, call for quote. Otherwise, provide a high quality color print of which is an accurate representation of the color expectations.
- **Actual Printed Dispenser:** 1 or 2 color, standard inks with camera-ready art provided: **\$35.00 (G).**
 - **Laser Engraved: \$145.00 (G).**
 - **Paper Insert: \$116.00 (G)** 1–4 colors. (DS330 only)
 - **Lenticular Insert: \$750.00 (G)** 1–4 colors. (DS330 only)
 - **Sheet Printing: \$145.00 (G)** for 1 or 2 color sheet imprint, standard or PMS inks with camera-ready art provided. All others call for quote.

Prototypes:

Refer to Speculative Samples on page 40.

Refills:

Unprinted, 3" x 3" x 50- or 100-sheet Post-it® Note Pad and Post-it® Flag Refills are available retail.

Restrictions:

- No serif fonts for engraving.
- No customization on Black base.
- No screens or gradient images for ink printed dispenser top.
- No fluorescent, special mix, or metallic ink on paper or lenticular insert.
- Paper and lenticular inserts not available for DS100. ■