

CLOSE OUT!

SPECIAL OFFER!

Valid Through December 31, 2009
or While Supplies Last!



Custom Printed
Flip Top Holder
with Pad

*Delivers
Advertising
Impact!*

*Folds flat for
easy mailing!*



Post-it® Flip Top Holder with Pad

- Folds flat for easy mailing.
- Each comes assembled with one (1) unprinted 3" x 4" 50-sheet Post-it® Note Pad.
- Imprint holder in up to 4 colors!

**Take \$1.00 OFF
every unit you buy!***

Post-it® Flip Top Holder with Pad — CLOSE OUT		Quantity Price Per Unit						
Item #	Description and Imprint	250-499	500-999	1000-2499	2500-4999	5000-9999	Set-up	
PCFTH34	Price includes: Holder and one (1) unprinted 3" x 4" 50-sheet pad, Standard white paper.							
	1-4 color imprint on holder	1.98 2.98	1.79 2.79	1.35 2.35	1.12 2.12	.72 1.72	BBCCD FREE	

Allow 8 days production time. Standard catalog upcharges apply.
Holder Imprint Area: 1-1/8" x 4". No bleeds on holder. No PMS or metallic inks. PMS colors will be converted to the closest possible 4-color match. Exact color match cannot be guaranteed. Note Pad printing available, call for production time and quote.
*Minimum order of 250 units.

Be sure to reference
offer #P4017
on purchase order

**8 Day
Production
Time!**

Exclusive offer valid through December 31, 2009 or while supplies last.