



Custom Printed  
Notes

# 24-Hour Service



No  
Additional  
Charge!

OUR MOST POPULAR SIZES!  
3" x 4" 25 or 50 sheet    4" x 6" 25 or 50 sheet

Post-it® Custom Printed Notes		Quantity Price per Unit							
	Item #	Sheet Count	Imprint	500	1000	2500		Units/Case Wt.	
<b>3" x 4"</b> Actual size: 2-7/8" x 4" Imprint area: 2-1/2" x 3-5/8" Adhesive side: 4" edge	PD34P-25	25-sheets	4-color process*	.69	.40	.36	bcc	640/21 lbs.	
	PD34P-50	50-sheets	4-color process*	.94	.71	.65	bcc	320/21 lbs.	
<b>4" x 6"</b> Actual size: 4" x 5-3/4" Imprint area: 3-5/8" x 5-3/8" Adhesive side: 4" edge	PD46P-25	25-sheets	4-color process*	1.02	.71	.65	bcc	320/21 lbs.	
	PD46P-50	50-sheets	4-color process*	1.60	1.22	1.05	bcc	160/21 lbs.	

Details for 24-Hour Service on other side.

\*For Full-Color printing, requests calling for a specific PMS color will be converted to the closest possible digital match. Exact color match cannot be guaranteed. No fluorescent, special mix or metallic inks.



No Additional Charge!



Custom Printed Lint Sheets Pocket Pack



Scotch® Custom Printed Lint Sheets Pocket Pack		Quantity Price per Unit**			
		Imprint	48	96	Units/Case Wt.
<b>PCLS4</b> Includes: 3" x 4" pad of 40 disposable lint remover sheets Imprint area on cover: 2-1/2" x 1-1/4" (Horizontal)	4-color process*	6.40	4.99	c	48/3.4 lbs. 96/6.4 lbs.
	High gloss finish (spot or full image) <b>add</b>	1.10	1.10	G	
	Copy change <b>add</b>	1.75	1.75	G	

U.S. Products and Pricing. Quantity price listed per unit.  
 \*For Full-Color printing, requests calling for a specific PMS color will be converted to the closest possible digital match. Exact color match cannot be guaranteed. No fluorescent, special mix or metallic inks.



Custom Printed Flag+ Writing Tools



Post-it® Custom Printed Flag+ Writing Tools		Quantity Price per Unit**			
		Imprint	40	80	Units/Case Wt.
<b>WIPFC4</b> Flag+ Pen <b>WIHFC4</b> Flag+ Highlighter Diameter: 5/8" Length: 5-1/2" Flags: (50) 3/8", unprinted Imprint area on barrel: 1-1/2" x 3/8" <b>White barrel only.</b>	4-color process*	4.29	2.99	c	40/1.8 lbs. 80/2.12 lbs.
	Copy change <b>add</b>	1.55	1.55	G	

		Imprint	100-249	250-499	500-999	1000-1999	2000-2499	Units/Case Wt.
<b>WIPHFC4</b> Flag+ Pen and Highlighter Combo Diameter: 1/2" Length: 5-11/16" Flags: (50) 5/16", unprinted Imprint area on barrel: 1-1/2" x 5/16"	4-color process*	3.38	2.69	2.59	2.29	2.19	BBBCC	500/19.25 lbs.
	Copy change <b>add</b>	quote	\$18.00/change			G		

24-Hour Service available by request only. Allow 24 hours after receipt of complete order or after approval of proof. Standard catalog upcharges apply. Available for PD34P-25, PD34P-50, PD46P-25, PD46P-50 catalog quantities only, and WIPFC4/WIHFC4 (White barrel only), PCLS4 low quantities only. 24-Hour Service must be in subject line of email to 24hour@3Mpromote.com. 24-Hour Service must be clearly stated on Purchase Order. 3M must receive complete Purchase Order and production-ready artwork by 1:00PM CST, your order ships the next business day. Call for units and case weights.

\*For Full-Color printing, requests calling for a specific PMS color will be converted to the closest possible 4-color process match. Exact color match cannot be guaranteed. No fluorescent, special mix or metallic inks.  
 \*\*Exact quantities only as listed above.

