



There's a note for that!™

most popular sizes

Post-it® Custom Printed Notes have MORE room to write and MORE room for your logo. Save on 3" x 4" and 4" x 6" size notes!



Post-it® Custom Printed Notes: **Exact Color Match** — Promotional Pricing Through December 31, 2009

When color matching to your logo or theme is critical, this program is ideal for you.

Item #	Size	Imprint	500	1000	2500	Set-Up
32	3" x 4" 25-sheets	1 spot color	.48	.39	.32	C FREE
		2 spot color	.69	.49	.38	
35	3" x 4" 50-sheets	1 spot color	.79	.70	.57	C FREE
		2 spot color	.99	.75	.62	
42	4" x 6" 25-sheets	1 spot color	.78	.68	.56	C FREE
		2 spot color	1.02	.78	.66	
45	4" x 6" 50-sheets	1 spot color	1.22	.99	.89	C FREE
		2 spot color	1.42	1.08	.94	
FREE PMS Ink Color Match*			FREE	FREE	FREE	FREE

Be sure to reference offer #P4076 on your P.O.

Allow 7 days production time. Quantity Price Per Pad. Wide array of standard paper colors available.

Add 15% upcharge for neon/ultra papers. Add 10% upcharge for recycled papers.

Special 5-Day Production Time on 500-2499 Quantity: Standard White paper only;


standard inks only. ●●●●●●●●

Items 32 and 35: **Actual size:** 2-15/16" x 4" **Imprint Area:** 2-9/16" x 3-5/8" **Adhesive side:** 4" edge

Items 42 and 45: **Actual size:** 4" x 5-7/8" **Imprint Area:** 3-5/8" x 5-1/2" **Adhesive side:** 4" edge

*No additional charge except fluorescent, special mix and metallic: \$285.00 (G) per color.

Choose a Paper Color:



Special Offers Valid September 1 – December 31, 2009

3M, Post-it, the color Canary Yellow and the statement "There's a note for that!" are trademarks of 3M. © 3M 2009. All Rights Reserved.

